

Collection Instructions for the BD Vacutainer Urine Collection Kit

1. Peel back protective sticker on top of the urine cup to expose rubber covered cannula.



2. Push C & S Preservative Tube (gray top tube) into the integrated transfer port.
 - Hold in position until flow stops.
 - Remove top.
 - Shake tube vigorously.



3. Push the UA Preservative Tube (red/ yellow top tube) into integrated transfer port.
 - Hold in position until flow stops.
 - Remove tube.
 - Invert UA Preservative Tube 8 – 10 times to mix the sample.



4. Place the protective sticker back over the integrated transfer port.
Label both filled tubes with the full patient name, date of birth, and the collection date/time.
5. Label the urine cup with the patient's full name and date of birth and submit urine cup and both tubes to the laboratory.